

Equine Standing Multidetector Computed Tomography of the Distal Thoracic Limb and Tarsus Has a Lower Cumulative Radiation Dose than Digital Radiography

Gaida et al. (2025), in Veterinary Radiology & Ultrasound

Products

Equine Standing Multidetector Computed Tomography of the Distal Thoracic Limb and Tarsus Has a Lower Cumulative Radiation Dose than Digital Radiography.

Hospital / Authors

Julia L. Gaida, Tim Steinberg, Susanne M. Stieger-Vanegas, Roswitha Merle, Christoph J. Lischer

Tierklinik Luesche GmbH, Bakum, Germany; Department of Clinical Sciences, Carlson College of Veterinary Medicine, Oregon State University, Corvallis, USA; Clinic for Horses, General Surgery and Radiology, School of Veterinary Medicine, Freie Universität Berlin, Berlin, Germany

Clinical Background

This study tests the hypothesis that MDCT radiation exposure to veterinary technicians will be equal to or lower than DR, challenging the belief that MDCT produces higher radiation exposure.

Aim of Study

To compare radiation exposure during MDCT and DR scans of the equine distal thoracic limb and tarsus, hypothesizing that MDCT will result in equal or lower exposure.

Cohort Study

24 Warmblood horses (12 each for MDCT and DR) were included, with imaging performed on the distal thoracic limb and tarsus for lameness or prepurchase exams between May 2021 and June 2022.

Results

- MDCT provided lower radiation exposure for most dosimeter locations (**thyroid**, **hands**, **feet**) compared to DR, supporting our hypothesis that radiation to veterinary staff during standing MDCT scans would be at least equal to, if not lower than, DR imaging of the same body area.
- The highest radiation exposure for MDCT occurred at the **gonads** during distal thoracic limb imaging, but overall exposure remained within safe limits.
- The average cumulative dose for MDCT was 0.61 μSv, while for DR, it was 1.43 μSv.
- MDCT scans took longer (approximately 20 minutes) compared to DR (around 3 minutes), which increased overall technician exposure, but this does not outweigh the lower cumulative dose.

Summary

- The findings support the hypothesis that MDCT is at least as safe, if not safer, than DR concerning radiation exposure to veterinary staff.
- The study emphasizes radiation protection strategies, such as increasing distance from the radiation source and using shielding techniques.
- MDCT, with its setup and protective measures, reduces radiation exposure to veterinary technicians compared to **DR**.
- Continued improvements in shielding and staff positioning are recommended to minimize radiation risks.
- MDCT offers a safer alternative for equine imaging, resulting in lower cumulative radiation doses compared to DR, ensuring better safety for veterinary staff.

<u>Link to paper</u>

Source: Gaida, J. L., Steinberg, T., Stieger-Vanegas, S. M., Merle, R., & Lischer, C. J. (2025). DOI: 10.1111/vru.13016, in Veterinary Radiology & Ultrasound.